USER INSTRUCTIONS

WIDEX BTE R D CHARGE n CLEAN

WPC201 model

GB
SC
DE
DK
ES
FR
IS
IT
NL
PL
PT


WIDEX®
HIGH DEFINITION HEARING

GB	USER INSTRUCTIONS	3
SC	用户使用手册	25
DE	BEDIENUNGSANLEITUNG	46
DK	BRUGSANVISNING	70
ES	INSTRUCCIONES DE USO	92
FR	MODE D'EMPLOI	116
IS	NOTKUNARLEIÐBEININGAR	140
IT	LIBRETTO D'ISTRUZIONI	162
NL	GEBRUIKSAANWIJZING	186
PL	INSTRUKCJA OBSŁUGI	210
PT	INSTRUÇÕES DE UTILIZAÇÃO	234
FCC	AND ISSED STATEMENTS	258

USER INSTRUCTIONS

IMPORTANT SAFETY INFORMATION.....	6
YOUR CHARGER.....	8
Welcome.....	8
The charger at a glance.....	8
Operating, storage and transport conditions.....	9
OPERATING THE CHARGER.....	10
Getting ready for use.....	10
Power supply.....	10
Power on.....	11
Charging the hearing aid.....	11
Charging indicators.....	12
After charging.....	14
Maintenance.....	14
Overview of indicator lights.....	16
Charging.....	16
Maintenance.....	17
Cleaning.....	18
TROUBLESHOOTING.....	19
REGULATORY INFORMATION.....	21
Intended use.....	21

EU directives.....	21
Directive 2014/53/EU.....	21
Information regarding disposal.....	22
SYMBOLS.....	23

IMPORTANT SAFETY INFORMATION

Read these pages carefully before you start using your charger.

Keep the charger and other small parts away from children, pets, and mentally disabled people.

Never expose the charger to temperatures above 60° C (140° F), to prolonged sunlight, open flames or ignition sources.

Place the charger upright on a flat, stable surface when charging, with adequate space around the charger. Do not cover or place near flammable substances. This also applies when using a power bank.

If cleaning is required, disconnect the charger from the power supply

Do not drop. Internal components may get damaged

Never try to open or repair the charger yourself. Contact your hearing care professional if you experience defects.

Do not use obviously damaged devices and return them to point of sale

Use only a WSAUD charger with your hearing aids. Use of a charger not recommended by WSAUD A/S may compromise the system or result in dangerous situations.

Never leave the charger in direct sunlight and do not immerse it in water.

Never insert other objects than your hearing aids in the charger.

Use only a WSAUD power adapter and cable with your charger.

Do not use when the status indicators are red

If despite the safety system, the ultraviolet (UV) light remains on when the charger lid is open, do not look directly at the UV light and contact your hearing care professional.

Do not use if charger is broken, and ultraviolet light is exposed.

YOUR CHARGER

Welcome

Thank you for choosing this product. Your new CHARGE n CLEAN device is intended to detect, charge and maintain your MBR3D hearing aid.


Read these instructions carefully before you start using this device.

NOTE

Your charger may not look exactly like the one illustrated in this booklet. We reserve the right to make any changes we consider necessary.

The charger at a glance

The charger consists of a charging unit for each hearing aid, a maintenance compartment for the receivers and a light indicator with three indicator lights that show the charging status of the hearing aids and the status of the maintenance cycle.

Operating, storage and transport conditions

Your charger is designed for indoor use. It is sensitive to extreme conditions such as high temperatures, and it should not be operated, stored or transported in direct sunlight.

Operate the charger within the following conditions:

	Minimum	Maximum
Temperature	0° C (32° F)	35° C (95° F)
Humidity	20% rH	75% rH

Store and transport the charger within the following conditions:

	Minimum	Maximum
Temperature	-20° C (-4° F)	60° C (140° F)
Humidity	5% rH	90% rH
Atmospheric pressure	750 mbar	1060 mbar

You can find technical data sheets and additional information on your charger on www.widex.com.

OPERATING THE CHARGER


Getting ready for use

Power supply

Power your charger using the included power supply. Use only the included cable and power supply.


Connect the cable to the charger:

Power on


When the charger is connected to the power supply, the indicator lights **flash white 5 times** and then turn off:

The charger is now in standby mode and ready to charge your hearing aid.

Charging the hearing aid

Gently place your hearing aids in their respective charging wells and the receivers in the maintenance compartment. The hearing aids turn off automatically. **Then close the lid.**

Left hearing aid

Right hearing aid

NOTE

The charger is not intended for hearing aids with anchors.


When connected to the power supply, the charger can charge two hearing aids simultaneously. Depending on the state of the hearing aids, each hearing aid can also be charged individually. If one hearing aid is fully charged and placed in the charger, it will not charge, but the other hearing aid, which is not fully charged, will be charged.

Charging your hearing aid fully will take approximately **4 hours**. If you charge your hearing aid for half an hour, you will have approximately 4 hours of use.

While charging the hearing aids, the maintenance (Ultraviolet-C treatment and drying) is carried out. The maintenance cycle takes 30 minutes.


Charging indicators

While the charger is detecting the hearing aid, the light indicator turns **solid yellow**:


Once the charger has detected the hearing aid, the charging will begin.


The light indicator will **flash green** to indicate that the hearing aid is low on battery and less than **40% charged**:

The light indicator will **glow green** to indicate that the hearing aid is more than **40% charged**:

The light indicator will **turn solid green** to indicate that the hearing aid is **100% charged**:

After charging


Avoid removing the hearing aids from the charging wells by pulling the earwire, as this could damage the connection between the earwire and the hearing aid. Instead, put your fingers around the hearing aid when removing it from the charger.

When you remove the hearing aid from the charger, the hearing aid turns on after a few seconds. Remember to turn it off if you are not going to use it right away.


Immediately after charging, the hearing aid may temporarily feel warm. If you are not comfortable with this, wait a few minutes for the hearing aid to reach room temperature, before putting the hearing aid on.

Maintenance

The light indicator will **flash white** to indicate that the maintenance cycle is in progress (Ultraviolet-C treatment and drying):

The light indicator turns **solid white** to indicate the end of the maintenance cycle:


If the lid is opened during the maintenance cycle, the cycle is paused but the charging continues. The maintenance cycle will restart when the lid is closed.

Maintenance pause
(**solid yellow** light)

Overview of indicator lights

Charging

Indicator lights	Status
Solid yellow 	The charger is detecting the hearing aid.
Flashing green 	The hearing aid is charging. Low battery, less than 40%.
Glowing green 	The hearing aid is charging. Battery more than 40%.
Solid green 	The hearing aid is fully charged. Battery 100%.
One indicator flashing red 	An error occurred during charging and the charging has stopped. See Troubleshooting.

Indicator lights

Status

Both indicators flashing red


An error occurred during charging and the charging has stopped.
See Troubleshooting.

Flashing yellow

The temperature of the hearing aid has become too high during charging and charging is paused.
See Troubleshooting.

Maintenance

Indicator lights

Status

Flashing white

The maintenance cycle is in progress (Ultraviolet-C treatment or drying).

Solid white

The maintenance cycle is complete.

Solid yellow

The maintenance cycle is paused but the charging cycle continues.

Indicator lights	Status
Flashing red 	An error occurred during the maintenance and the cycle has stopped. See Troubleshooting.
Flashing yellow 	The maintenance cycle is paused due to overheating (fan still running). See Troubleshooting.

Cleaning

The charger does not generally require cleaning. If necessary, carefully wipe off the charger with a soft dry cloth.

If you notice any dust in the charging well, gently blow it clean. If this is not enough, gently clean the charging well with a cotton bud.

To remove earwax from the charger, carefully wipe with a soft cloth with a bit of mild household detergent.

Never use sharp/pointy objects to clean the charger.

TROUBLESHOOTING

These pages contain advice on what to do if your charger does not work as expected.

Light indicators	Possible cause	Solution
No lights	The charger is not connected.	Make sure that the USB cable is correctly connected to the charger and the wall socket using the included power adaptor.
	The hearing aid is not placed correctly in the charger.	Place the hearing aid correctly in the charger.
Right and/or left indicators flashing yellow 	The temperature of your surroundings or of the hearing aid exceeds a critical limit and the charging is paused.	The charging will continue when the temperature of your surroundings or the hearing aid has dropped. If this happens continually, find a cooler place to charge your hearing aid.
Right or left indicator flashing red 	The charger cannot recognise the hearing aid.	Remove the incompatible hearing aid.

Light indicators	Possible cause	Solution
	The hearing aid has a critical error.	Remove the hearing aid, disconnect from the power for 10 seconds and then try again. If the charger is still flashing red, contact your hearing care professional.
Central indicator flashing red 	The maintenance cycle has a critical error.	Remove the hearing aid, disconnect from the power for 10 seconds and then try again. If the charger is still flashing red, contact your hearing care professional.
Right and left indicators flashing red 	The charger has a critical error.	Remove the hearing aid, disconnect from the power for 10 seconds and then try again. If the charger is still flashing red, contact your hearing care professional.

REGULATORY INFORMATION

Intended use

The WIDEX charger WPC201 model is to be used for charging the MBR3D hearing aid.

EU directives

Directive 2014/53/EU

Hereby, WSAUD A/S declares that the WPC201 is in compliance with the essential requirements and other relevant provisions of Directive 2014/53/EU.

The WPC201 contains a radio transmitter operating at:
6.78 MHz, 13.4 dB μ A/m @ 3 m.

A copy of the Declaration of Conformity according to 2014/53/EU can be found at

www.widex.com/doc

N26346

Information regarding disposal

Do not dispose of hearing aids, hearing aid accessories and chargers with ordinary household waste.

Hearing aids, hearing aid accessories and chargers should be disposed of at sites intended for electrical waste and electronic equipment, or given to your hearing care professional for safe disposal.

For returning or shipping the devices, please contact your hearing care professional.

Shipment of charger and hearing instruments should follow to local regulations.

SYMBOLS

Symbols commonly used by WIDEX in medical device labelling (labels/IFU/ etc.)

Symbol	Title/Description
	Manufacturer The product is produced by the manufacturer whose name and address are stated next to the symbol. If appropriate, the date of manufacture may also be stated.
	Catalogue number The product's catalogue (item) number.
	Consult instructions for use The user instructions contain important cautionary information (warnings/precautions) and must be read before using the product.
	Warning Text marked with a warning symbol must be read before using the product.
	WEEE mark “Not for general waste”. When a product is to be discarded, it must be sent to a designated collection point for recycling and recovering to prevent the risk of harm to the environment or human health as a result of the presence of hazardous substances.

Symbol	Title/Description
	CE mark The product is in conformity with the requirements set out in European CE marking directives.
	RCM mark The product complies with electrical safety, EMC and radio spectrum regulatory requirements for products supplied to the Australian or New Zealand markets.
	FCC logo The product is compliant with relevant USA 47 CFR Rules. FCC – Federal Communications Commission.
	Medical device Indication that the device is a medical device.

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FCC AND ISED STATEMENTS

Federal Communications Commission Statement

This device complies with part 15/part 18 of the FCC Rules. Operation is subject to the following two conditions:

- (1) This device may not cause harmful interference, and
- (2) this device must accept any interference received, including interference that may cause undesired operation.

NOTE:

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.

- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

NOTE:

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

Changes or modifications to the equipment not expressly approved by WSAUD A/S could void the user's authority to operate the equipment.

ISED Statement / Déclaration d'ISED

CAN RSS-216/CNR-216

CAN ICES-1/NMB-1

This device contains licence-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada's licence-exempt RSS(s). Operation is subject to the following two conditions:

- (1) This device may not cause interference.
- (2) This device must accept any interference, including interference that may cause undesired operation of the device.

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ISED RADIATION EXPOSURE STATEMENT:

This equipment complies with ISED RF radiation exposure limits set forth for an uncontrolled environment. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

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